

Page: 1

Date: 01/19/22 at 10:03 AM

## Aqua Solutions, Inc. 6913 Highway 225 Deer Park, TX, 77536 281/479-2569

## **Certificate of Analysis / QC Results**

Customer PO #:

Packaged Product: 6520			
Test	Target/IOM		Result
PREPARED TO FORMULATION	Target/UOM ON YES	Range	YES
FILE	YES		120
Sample ID	120		
	2		
Tested 01/19/2	2		
Tester RIRO	•		
Recorded 01/19/2			0.50
STIR FOR A NUMBER OF HOUR		2.00 - 3.50	3.50
	YES		
Sample ID			
Tested 01/19/2	2		
Tester RIRO			
Recorded 01/19/2	2		
APPEARANCE AND COLOR	CLEAR WATER WHITE		PASS
	LIQUID		
	APPEARANCE		
Sample ID			
Tested 01/19/2	2		
Tester RIRO			
Recorded 01/19/2	2		
INSTRUMENTS USED DURING	INSTRUMENT		T-D3,M55,VF3-J11
PREPARATION			. 20,00, 0 0
Sample ID			
Tested 01/19/2	2		
	2		
Tester RIRO	2		
Recorded 01/19/2	0.82686	0.82547 - 0.82824	0.82719
Actual Density		0.62547 - 0.62624	0.82719
	GM/ML		
Sample ID			
Tested 01/19/2	2		
Tester RIRO			
Recorded 01/19/2			
VERIFICATION	PASS		PASS
	PASS		
Sample ID			
Tested 01/19/2	2		
Tester RIRO			
Recorded 01/19/2	2		
EXPIRATION DATE	MM/DD/YY		01/30/24
	MM/DD/YYYY		
Sample ID			
Tested 01/19/2	2		
Tester RIRO	_		
Recorded 01/19/2	2		
SHELF LIFE OF PRODUCT	SHELF LIFE		1 YEAR
	TIME		LILAIT
Comple ID	i iiviL		
Sample ID	2		
Tested 01/19/2	2		
Tester RIRO	2		
Recorded 01/19/2	₹		

Page: 2

Date: 01/19/22 at 10:03 AM

## Aqua Solutions, Inc. 6913 Highway 225 Deer Park, TX, 77536 281/479-2569 Certificate of Analysis / QC Results

Customer PO #:

Lot # 2011435 Made 01/19/22

This is to certify that the product listed above has been prepared according to the agreed upon formulation. The solutions producer has a Quality Management System which governs each step of the manufacturing process to insure the production of a consistent product. Traceability from the producer's lot numbers to the original manufacturer's lot number is maintained. Certificates of Analysis for the raw materials used to prepare this product are available upon request.

The weights and/or volumes used to prepare this product are N.I.S.T. Traceable. All balances used in the preparation of product are calibrated daily against N.I.S.T. Traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is N.I.S.T. Traceable and certified as meeting Class A specifications.

Unless otherwise agreed to, all chemicals used in the preparation of this product are Reagent ACS grade.

**Quality Manager**