

Recorded 07/18/22

Page: 1

Date: 08/13/22 at 8:47 AM

Aqua Solutions, Inc. 6913 Highway 225 Deer Park, TX, 77536 281/479-2569 Certificate of Analysis / QC Results

Customer PO #:

Packaged Product:				
	EDTA TI	TRANT 0.05 MOLAR (0.11)	
Test		Target/UOM	Range	Result
PREPARED TO FORMU	LATION ON	YES		YES
FILE		YES		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
Actual Molarity		0.0500	0.0495 - 0.0505	0.0496
		MOL/L		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
QCSP4:STANDARDIZAT	ION EDTA	YES		YES
SOLUTIONS		QCSP4		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
N.I.S.T. Traceable to SR	M 915	YES		YES
		YES		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
APPEARANCE AND COLOR		CLEAR LIQUID		PASS
		APPEARANCE		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
INSTRUMENTS USED DURING		INSTRUMENT		M103,M55,(O)T-141,BID03828
PREPARATION				PIDM81
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
EXPIRATION DATE		MM/DD/YY		07/30/24
		MM/DD/YYYY		
Sample ID				
Tested	07/18/22			
Tester	CCOA			
Recorded	07/18/22			
Actual Density		1.00749	1.00710 - 1.00789	1.00730
		GM/ML		
Sample ID				
Tested	07/18/22			
Tester	CCOA			

Page: 2

Date: 08/13/22 at 8:47 AM

Aqua Solutions, Inc. 6913 Highway 225 Deer Park, TX, 77536 281/479-2569 Certificate of Analysis / QC Results

Customer PO #:

SHELF LIFE OF PRODUCT 2 YEAR 2 YEAR TIME

Sample ID

Tested 07/18/22
Tester CCOA
Recorded 07/18/22

Lot # 2071242 Made 07/18/22

This is to certify that the product listed above has been prepared according to the agreed upon formulation. The solutions producer has a Quality Management System which governs each step of the manufacturing process to insure the production of a consistent product. Traceability from the producer's lot numbers to the original manufacturer's lot number is maintained. Certificates of Analysis for the raw materials used to prepare this product are available upon request.

The weights and/or volumes used to prepare this product are N.I.S.T. Traceable. All balances used in the preparation of product are calibrated daily against N.I.S.T. Traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is N.I.S.T. Traceable and certified as meeting Class A specifications.

Unless otherwise agreed to, all chemicals used in the preparation of this product are Reagent ACS grade.

Quality Manager